

Installation Instructions on Medonic M51

A detailed description of reagent installation / reagent change can be found in the User Manual of the analyzer.



Intended Use

The Boule Medonic M51-D Diluent, Medonic M51-L1 Lyse, and Medonic M51-L2 Lyse reagents are intended for use with the Medonic M51 hematology analyzer.



Storage and Stability

The Boule Medonic M51 reagents are stable at least up until the expiration date as shown on labels, when stored at 2–30 °C. Open container stability is 60 days if storage recommendations are followed.

Instructions

For detailed instructions, please refer to the User Manual of the analyzer.

- For a proper performance of the analyzer, DO NOT mix Boule Medonic M51 reagents with non-Boule products.
- It is recommended that the Boule Medonic M51 reagents are placed at the same level as the analyzer or below.
- When a reagent is running out, an error message is prompted on the screen. Perform the following actions:
 - a) Select the error message and then select the error in the popup dialog box and click “Remove Error”.
 - b) A dialog box pops up. “Put the RFID card on the card reader”.
 - c) Place the RFID card (found on the reagent container) on the RFID symbol on the front of the analyzer. The card reader beeps and a popup dialog box showing successful reagent settings appears. Click OK in the popup dialog box.
 - d) For the Medonic M51 Lyse reagent, unscrew the cap of the reagent container, put the aspiration tube into the container and screw the cap. Note the color coding (the Medonic M51-L1 Lyse is green and the Medonic M51-L2 Lyse is blue).
For the Medonic M51-D Diluent reagent, put the float sensor into the container and screw the cap. The float sensor should be kept as vertical as possible during installation.
 - e) Run a background check and ensure that the results are within the allowed limits.

Notes

- Prevent contamination of the reagents with dust, etc. by use of the cap (delivered with the analyzer).
- The Diluent container can be reused as a waste container.
- The Diluent cardboard box can be recycled after removal of the LDPE container. Emptied LDPE containers and HDPE Lyse containers can be rinsed with water and recycled.

Ordering information and service

For international customers contact your local distributor for orders and support.

REF	Description	Packaging
1504510	Medonic M51-D Diluent	20 L
1504511	Medonic M51-L1 Lyse	200 mL
1504512	Medonic M51-L2 Lyse	500 mL