

# Boule Cal

CAL

IVD



Calibrator

Lot	
22405-34	2024-08-19

Open vial stability 5 days <sup>1</sup>

## Medonic

Parameter <sup>2</sup>	RBC 10 <sup>12</sup> /L	MCV fL	Plt 10 <sup>9</sup> /L	MPV fL	WBC 10 <sup>9</sup> /L	Hgb g/dL
Medonic M-series / M32	4.21	81.4	238	11.5	8.6	12.2
Range ± <sup>3</sup>	0.10	2.0	15	1.5	0.3	0.2

## Swelab

Parameter <sup>2</sup>	RBC 10 <sup>12</sup> /L	MCV fL	Plt 10 <sup>9</sup> /L	MPV fL	WBC 10 <sup>9</sup> /L	Hgb g/dL
Swelab Alfa / Alfa Plus	4.21	82.4	243	12.4	8.6	12.2
Range ± <sup>3</sup>	0.10	2.0	15	1.5	0.3	0.2

	1	2	3
bg	Трайност на отворен флакон	5 дни	Параметър
hr	Stabilnost otvorene bočice	5 dani	Parametar
cs	Stabilita otevřené ampulky	5 dní	Parametr
da	Åben rør stabilitet	5 dage	Parameter
et	Stabiilsus avatud vialis	5 päevade arv	Parameeter
fr	Stabilité en flacon ouvert	5 jour	Paramètre
de	Stabilität geöffneter Flaschen	5 tage	Parameter
el	Διατήρηση δείγματος μετά την αποσφράγιση	5 ημέρα	Παράμετρος
it	Stabilità della fiala aperta	5 giorni	Parametro
lv	Atvērt flakonu stabilitāti	5 dienas	Parametrs
lt	Stabilumas atidarius buteliuka	5 d.	Parametras
no	Åpen rør/glass stabilitet	5 dager	Parameter
pl	Trwałość po otwarciu fioleki	5 Liczba dni	Kryterium
pt	Estabilidade após abertura do frasco	5 dias	Parâmetro
ro	Termenul de valabilitate al fiolei desfăcută	5 zile	Parametru
ru	Стабильность вскрытого флакона	5 дней	Параметр
sr	Proizvod stabilan nakon otvaranja pakovanja	5 dana	Parametar
sk	Stabilita otvorenej injekčnej liekovky	5 dní	Parameter
sl	Stabilnost odprte stekleničke	5 dni	Parameter
es	Estabilidad de la cápsula abierta	5 días	Parámetro
sv	Hållbarhet för öppnad flaska	5 dagar	Parameter
tr	Açık şişe stabilitesi	5 gün	Parametre
			Обхват
			Raspon
			Rozsah
			Område
			Vahemik
			Intervalle
			Bereich
			Αναμενόμενο εύρος α
			Intervallo
			Diapazons
			Intervalas
			Område
			Zakres
			Intervalo
			Interval
			диапазон
			Opseg
			Rozsah
			Območje
			Intervalo
			Intervall
			Aralığı

[www.boule.com](http://www.boule.com)

Boule Medical AB  
Domnarvsgatan 4  
SE-163 53 Spånga, Sweden

25123-6

1504055

## For Medonic M-series / M32 with barcode reader

### M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

### M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series

Lot: 22405-34

Prod.: Cal

ID: 2240534+

Exp.: 19-Aug-2024



## For Swelab Alfa / Alfa Plus with barcode reader

### Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

### Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa  
Lot: 22405-34  
Prod.: Cal  
ID: 2240534+  
Exp.: 19-Aug-2024

