

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22405-1K	2024-10-04
Low	22405-11	2024-10-04
Normal	22405-12	2024-10-04
High	22405-13	2024-10-08

MEDONIC

CONTROL		L	
LOT		22405-11	
Medonic M-series/M32	Mean ²	± Range ³	
RBC	10 ¹² /L	2.30	0.12
MCV	fL	67.3	5.0
RDW	%	15.1	4.0
Hct	%	15.5	2.5
Plt	10 ⁹ /L	85	15
MPV	fL	11.4	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.7	0.3
MCH	pg	24.8	1.5
MCHC	g/dL	36.8	3.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	44.3	7.0
Mid	%	5.9	5.0
Gran	%	49.8	7.0

CONTROL		N	
LOT		22405-12	
	Mean ²	± Range ³	
	4.20	0.18	
	86.0	5.0	
	13.8	4.0	
	36.1	3.0	
	236	30	
	11.6	1.5	
	8.7	0.6	
	12.6	0.4	
	30.0	1.5	
	34.9	3.0	
	3.9	0.6	
	0.6	0.5	
	4.2	1.0	
	44.8	5.0	
	7.3	5.0	
	47.9	5.0	

CONTROL		H	
LOT		22405-13	
	Mean ²	± Range ³	
	5.15	0.22	
	94.8	5.0	
	14.2	4.0	
	48.8	3.5	
	525	60	
	12.0	1.5	
	20.9	1.8	
	16.5	0.5	
	32.0	1.5	
	33.8	3.0	
	12.6	1.5	
	1.2	1.0	
	7.1	2.5	
	60.3	5.0	
	5.6	5.0	
	34.1	5.0	

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dni (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αεροσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonū stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelværdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
Lot: 22405-11
Prod.: Con-Diff Low
ID: 2240511+
Exp.: 4-Oct-2024

Instr.: M-series
Lot: 22405-12
Prod.: Con-Diff Norm
ID: 2240512+
Exp.: 4-Oct-2024

Instr.: M-series
Lot: 22405-13
Prod.: Con-Diff High
ID: 2240513+
Exp.: 8-Oct-2024



Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22405-1K	2024-10-04
Low	22405-11	2024-10-04
Normal	22405-12	2024-10-04
High	22405-13	2024-10-08

SWELAB		CONTROL	L
Swelab Alfa/Alfa Plus		LOT	22405-11
		Mean	± Range
RBC	10 ¹² /L	2.31	0.12
MCV	fL	68.5	5.0
RDW	%	14.4	4.0
Hct	%	15.8	2.5
Plt	10 ⁹ /L	86	15
MPV	fL	12.5	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.7	0.3
MCH	pg	24.7	1.5
MCHC	g/dL	36.1	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.6	0.8
Lymph	%	47.7	7.0
Mid	%	4.2	5.0
Gran	%	48.1	7.0

CONTROL	N
LOT	22405-12
Mean	± Range
4.22	0.18
87.5	5.0
12.5	4.0
36.9	3.0
237	30
12.6	1.5
8.6	0.6
12.6	0.4
29.9	1.5
34.1	3.0
4.2	0.6
0.4	0.5
4.0	1.0
48.6	5.0
4.6	5.0
46.8	5.0

CONTROL	H
LOT	22405-13
Mean	± Range
5.13	0.22
96.3	5.0
13.3	4.0
49.4	3.5
538	60
13.8	1.5
20.5	1.8
16.3	0.5
31.8	1.5
33.0	3.0
13.1	1.5
0.7	1.0
6.7	2.5
63.8	5.0
3.4	5.0
32.8	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dni (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud viaalis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αλοσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonū stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelværdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 22405-11
Prod.: Con-Diff Low
ID: 2240511+
Exp.: 4-Oct-2024

Instr.: Swelab Alfa
Lot: 22405-12
Prod.: Con-Diff Norm
ID: 2240512+
Exp.: 4-Oct-2024

Instr.: Swelab Alfa
Lot: 22405-13
Prod.: Con-Diff High
ID: 2240513+
Exp.: 8-Oct-2024

