

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹



Control	Lot	
Tri-Level	32409-5K	2025-01-30
Low	32409-51	2025-01-30
Normal	32409-52	2025-01-31
High	32409-53	2025-02-05

MEDONIC

CONTROL	L
LOT	32409-51

Medonic M-series/M32		Mean ²	± Range ³
RBC	10 ¹² /L	2.24	0.12
MCV	fL	66.3	7.0
RDW	%	15.6	4.5
Hct	%	14.9	3.0
Plt	10 ⁹ /L	80	15
MPV	fL	11.1	1.5
WBC	10 ⁹ /L	3.5	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.6	2.0
MCHC	g/dL	36.9	4.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	45.7	7.0
Mid	%	5.5	5.0
Gran	%	48.8	7.0

CONTROL	N
LOT	32409-52

		Mean ²	± Range ³
		4.20	0.18
		81.7	7.0
		14.8	4.5
		34.3	3.5
		231	30
		11.3	1.5
		8.7	0.6
		12.0	0.4
		28.6	2.0
		35.0	4.0
		4.0	0.6
		0.6	0.5
		4.1	1.0
		46.5	5.0
		6.4	5.0
		47.1	5.0

CONTROL	H
LOT	32409-53

		Mean ²	± Range ³
		5.07	0.22
		90.5	7.0
		14.0	4.5
		45.9	4.0
		523	60
		11.6	1.5
		20.9	1.8
		15.8	0.5
		31.2	2.0
		34.4	4.0
		12.8	1.5
		1.2	1.0
		6.9	2.5
		61.2	5.0
		5.8	5.0
		33.0	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevrené ampulky 14 dny (cs)
 - da Åben r r stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 p evade (et)
 - fr Stabilit  en flacon ouvert 14 jour (fr)
 - de Stabilit t ge ffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημε ρα (el)
 - it Stabilit  della fiala aperta 14 giorni (it)
 - lv Atv rt flakonu stabilit ti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuk  14 dienas (lt)
 - no  pen r r/glass stabilitet 14 dager (no)
 - pl Trwa o c po otwarciu fioleki 14 Liczba dni (pl)
 - pt Estabilidade ap s abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desf cute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injek nej liekovky 14 dn  (sk)
 - sl Stabilita otvorenej injek nej liekovky 14 dni (sl)
 - es Estabilidad de la c psula abierta 14 dias (es)
 - sv H llbarhet f r  ppnad flaska 14 dagar (sv)
 - tr A ık ŐiŐe stabilitesi 14 g n (tr)

- 2
- Sredno
 - Sredstvo
 - Pr m r
 - Middelv rdi
 - Keskm ne
 - Moyenne
 - Mittelwert
 - M Ő 
 - Media
 - Noz m 
 - Vidurkis
 - Middelverdi
 - Őrednia
 - M dia
 - Medie
 - Srednee
 - Srednja vrednost
 - Priemer
 - Povpre no
 - Media
 - Medelv rde
 - Ortalama

- 3
- Ob vat
 - Raspon
 - Rozsah
 - Omr de
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμεν μενο  ρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Omr de
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Obmo je
 - Intervalo
 - Intervall
 - Aralıđı

25122-4



[http://www.medonic.se/
support/download](http://www.medonic.se/support/download)

Boule
Boule Medical AB
Domnarvsgatan 4
SE-163 53 Sp nga,
Sweden

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
 Lot: 32409-51
 Prod.: Con-Diff Low
 ID: 3240951+
 Exp.: 30-Jan-2025

Instr.: M-series
 Lot: 32409-52
 Prod.: Con-Diff Norm
 ID: 3240952+
 Exp.: 31-Jan-2025

Instr.: M-series
 Lot: 32409-53
 Prod.: Con-Diff High
 ID: 3240953+
 Exp.: 5-Feb-2025



Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹



Control	Lot	
Tri-Level	32409-5K	2025-01-30
Low	32409-51	2025-01-30
Normal	32409-52	2025-01-31
High	32409-53	2025-02-05

SWELAB

		CONTROL	L
		LOT	32409-51
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.27	0.12
MCV	fL	68.0	7.0
RDW	%	15.1	4.5
Hct	%	15.4	3.0
Plt	10 ⁹ /L	80	15
MPV	fL	12.3	1.5
WBC	10 ⁹ /L	3.5	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	24.7	2.0
MCHC	g/dL	36.4	4.0
Lymph#	10 ⁹ /L	1.7	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	48.4	7.0
Mid	%	4.0	5.0
Gran	%	47.6	7.0

		CONTROL	N
		LOT	32409-52
		Mean ²	± Range ³
		4.21	0.18
		83.1	7.0
		13.6	4.5
		35.0	3.5
		237	30
		12.5	1.5
		8.6	0.6
		12.1	0.4
		28.7	2.0
		34.6	4.0
		4.3	0.6
		0.4	0.5
		3.9	1.0
		50.5	5.0
		4.2	5.0
		45.3	5.0

		CONTROL	H
		LOT	32409-53
		Mean ²	± Range ³
		5.07	0.22
		92.0	7.0
		12.9	4.5
		46.6	4.0
		541	60
		12.8	1.5
		20.8	1.8
		15.8	0.5
		31.2	2.0
		33.9	4.0
		13.3	1.5
		0.7	1.0
		6.8	2.5
		63.8	5.0
		3.6	5.0
		32.6	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelværdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Orseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı

25122-4



www.swelab.com/
support/download



For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 32409-51
Prod.: Con-Diff Low
ID: 3240951+
Exp.: 30-Jan-2025

Instr.: Swelab Alfa
Lot: 32409-52
Prod.: Con-Diff Norm
ID: 3240952+
Exp.: 31-Jan-2025

Instr.: Swelab Alfa
Lot: 32409-53
Prod.: Con-Diff High
ID: 3240953+
Exp.: 5-Feb-2025

