

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	--	--
Low	22503-41	2025-08-04
Normal	22503-42	2025-08-04
High	22503-43	2025-08-06

MEDONIC

Medonic M-series/M32		CONTROL	L
		LOT	22503-41
		Mean ²	± Range ³
RBC	10 ¹² /L	2.28	0.12
MCV	fL	65.7	5.0
RDW	%	15.9	4.0
Hct	%	15.0	2.5
Plt	10 ⁹ /L	83	15
MPV	fL	11.1	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.1	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	43.1	7.0
Mid	%	6.0	5.0
Gran	%	50.9	7.0

CONTROL	N	LOT	22503-42	Mean ²	± Range ³
				4.23	0.18
				82.8	5.0
				13.4	4.0
				35.0	3.0
				240	30
				11.4	1.5
				8.3	0.6
				12.3	0.4
				29.1	1.5
				35.1	3.0
				3.9	0.6
				0.6	0.5
				3.8	1.0
				46.8	5.0
				7.2	5.0
				46.0	5.0

CONTROL	H	LOT	22503-43	Mean ²	± Range ³
				5.16	0.22
				93.1	5.0
				13.6	4.0
				48.0	3.5
				541	60
				11.6	1.5
				20.3	1.8
				16.5	0.5
				32.0	1.5
				34.4	3.0
				12.5	1.5
				1.2	1.0
				6.6	2.5
				61.6	5.0
				5.9	5.0
				32.5	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelværdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series	Instr.: M-series	Instr.: M-series
Lot: 22503-41	Lot: 22503-42	Lot: 22503-43
Prod.: Con-Diff Low	Prod.: Con-Diff Norm	Prod.: Con-Diff High
ID: 2250341+	ID: 2250342+	ID: 2250343+
Exp.: 4-Aug-2025	Exp.: 4-Aug-2025	Exp.: 6-Aug-2025

1:  B o u l e C B A V 2	 B o u l e C B A V 2	 B o u l e C B A V 2
2:  2250341+	 2250342+	 2250343+
3:  25080400001153030000	 25080400002153030000	 25080600003153030000
4:  10004030700000000000	 10004030200000000000	 10004030100000000000
5:  22812657500550315025	 42318828501230435030	 51622931501650548035
6:  08315111152411536730	 24030114152911535130	 54160116153201534430
7:  03403015040020301708	 08306039060060503810	 20318125150121006625
8:  15940431700605050970	 13440468500725046050	 13640616500595032550
9:  11110330410241854581	 11117130521186319659	 11117480289093787171

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	--	--
Low	22503-41	2025-08-04
Normal	22503-42	2025-08-04
High	22503-43	2025-08-06

Open vial stability 14 days ¹

SWELAB

Swelab Alfa/Alfa Plus	CONTROL		L
	LOT	22503-41	22503-41
	Mean ²	± Range ³	
RBC	10 ¹² /L	2.31	0.12
MCV	fL	67.0	5.0
RDW	%	15.1	4.0
Hct	%	15.5	2.5
Plt	10 ⁹ /L	86	15
MPV	fL	11.8	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.7	0.3
MCH	pg	24.7	1.5
MCHC	g/dL	36.8	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	47.1	7.0
Mid	%	3.3	5.0
Gran	%	49.6	7.0

CONTROL	N	
	LOT	22503-42
	Mean ²	± Range ³
	4.29	0.18
	84.1	5.0
	12.4	4.0
	36.1	3.0
	249	30
	12.0	1.5
	8.3	0.6
	12.4	0.4
	28.9	1.5
	34.3	3.0
	4.3	0.6
	0.3	0.5
	3.7	1.0
	50.9	5.0
	4.1	5.0
	45.0	5.0

CONTROL	H	
	LOT	22503-43
	Mean ²	± Range ³
	5.21	0.22
	94.5	5.0
	12.6	4.0
	49.2	3.5
	569	60
	12.3	1.5
	20.2	1.8
	16.4	0.5
	31.5	1.5
	33.3	3.0
	13.1	1.5
	0.7	1.0
	6.4	2.5
	64.8	5.0
	3.5	5.0
	31.7	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilnost otvorene bočice 14 dani (hr)
cs	Stabilita otevřené ampulky 14 dní (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud vialis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiołki 14 Liczba dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvorenej injekčnej liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 días (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

	Средно
	Sredstvo
	Průměr
	Middelverdi
	Keskmine
	Moyenne
	Mittelwert
	Μέση
	Media
	Nozīmē
	Vidurkis
	Middelverdi
	Średnia
	Média
	Medie
	Среднее
	Srednja vrednost
	Priemer
	Povprečno
	Media
	Medelvärde
	Ortalama

3

	Обхват
	Raspon
	Rozsah
	Område
	Vahemik
	Intervalle
	Bereich
	Αναμενόμενο εύρος α
	Intervallo
	Diapazons
	Intervalas
	Område
	Zakres
	Intervalo
	Interval
	диапазон
	Orseg
	Rozsah
	Območje
	Intervalo
	Intervall
	Aralığı

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
 Lot: 22503-41
 Prod.: Con-Diff Low
 ID: 2250341+
 Exp.: 4-Aug-2025

Instr.: Swelab Alfa
 Lot: 22503-42
 Prod.: Con-Diff Norm
 ID: 2250342+
 Exp.: 4-Aug-2025

Instr.: Swelab Alfa
 Lot: 22503-43
 Prod.: Con-Diff High
 ID: 2250343+
 Exp.: 6-Aug-2025

1:	 Boul e CBAV3	 Boul e CBAV3	 Boul e CBAV3
2:	 2250341+	 2250342+	 2250343+
3:	 25080400001153030000	 25080400002153030000	 25080600003153030000
4:	 20014018000000000000	 20014018000000000000	 20014018000000000000
5:	 23112670500570315525	 42918841501240436130	 52122945501640549235
6:	 08615118152471536830	 24930120152891534330	 56960123153151533330
7:	 03403016040010301708	 08306043060030503710	 20218131150071006425
8:	 15140471700335049670	 12440509500415045050	 12640648500355031750
9:	 11119614517492145075	 11118155567768479017	 11113401952262691941