

## Installation Instructions on M-series

A detailed description of reagent installation / reagent change can be found in the instrument User's Manual under section 2 Installation



## Intended use

The Boule Medonic M-series Diluent and Lyse reagents are intended for use with the Medonic M-series hematology analyzer systems

IVD

## Storage and stability

The Boule Medonic M-series reagents are stable at least up until the expiration date as shown on labels when stored at 4 to 35 °C. Open container stability 3 months if storage recommendations are followed.

## Instructions

1. For proper performance **DO NOT** mix Boule reagents with non-Boule products.
2. It is recommended that both the Diluent and the Lyse reagents are placed at the instrument level or below.
3. The interlocked reagent system displays indicator and warning messages to alert the operator when reagents are running low and need to be changed. When this occurs perform the following:
  - a. Select [MENU] to access the Main menu and then select [REAGENT SETUP].
  - b. Select [ENTER NEW REAGENT].
  - c. Scan Barcode 1 and then Barcode 2 on the reagent container. Press and hold the ON button on the barcode reader each time a barcode is scanned.
  - d. When all barcodes are entered a screen will display that reagent barcodes have been accepted.
  - e. Select [EXIT] to return to the Main menu.
  - f. Remove the cap and seal on the new reagent container.
  - g. Transfer the reagent level sensor from the used container to the new reagent container.
  - h. The analyzer is now ready to resume operation or analyze samples. No priming or fill cycle is necessary when putting on a new reagent container, if indicator and warning messages are followed.
  - i. Run three (3) BACKGROUND CYCLES to complete the reagent exchange

You can always refer to your Medonic M-series User's manual for reagent priming or changing.

**IMPORTANT!** Remember to verify calibration prior to running patient samples e.g. using Boule control material.

## Notes

- Prevent contamination of the reagents with dust, etc by use of the protective cap (delivered with the instrument).
- The cardboard box can be recycled after removal of the LPDE cubitainers. The LPDE container can be recycled after removal of any remaining reagent and a rinse with water.

## Ordering information and service

For international customers contact your local distributor for orders and support.

REF	Description	Packaging
1504122	Medonic M-series Diluent	20 L
1504123	Medonic M-series Lyse	5 L
1504128	Medonic M-series DualPack	1 L Lyse / 5 L Diluent
1504460	Medonic M-series Diluent, RFID	20 L
1504461	Medonic M-series Lyse, RFID	5 L
1504465	Medonic M-serie DualPack, RFID	1 L Lyse / 5 L Diluent



### Installation Instructions on M-series M32

A detailed description of reagent installation / reagent change can be found in the instrument User's Manual under section 2 Installation



### Intended use

The Boule Medonic M-series Diluent and Lyse reagents are intended for use with the Medonic M-series hematology analyzer systems

IVD

### Storage and stability

The Boule Medonic M-series reagents are stable at least up until the expiration date as shown on labels when stored at 4 to 35 °C. Open container stability 3 months if storage recommendations are followed.

### Instructions

1. For proper performance **DO NOT** mix Boule reagents with non-Boule products.
2. It is recommended that both the Diluent and the Lyse reagents are placed at the instrument level or below.
3. The interlocked reagent system displays indicator and warning messages to alert the operator when reagents are running low and need to be changed. When this occurs perform the following:
  - a. Select **Quick Functions Menu** and then select **Add Reagent**.
  - b. Scan in barcodes/Rfid tags on reagent box, and when all reagent barcodes/Rfid tags are entered a screen will display that reagent barcodes/Rfid tags have been accepted.
  - c. Select **Exit** to return to the **Quick Functions Menu**.
  - d. Remove the cap and seal on the new reagent container.
  - e. Transfer the reagent tube assembly from the used container to the new reagent container.
  - f. The analyzer is now ready to resume operation or analyze samples. No priming or fill cycle is necessary when putting on a new reagent container, if indicator and warning messages are followed.

You can always refer to your Medonic M-series M32 User's manual for reagent priming or changing.

**IMPORTANT!** Remember to verify calibration prior to running patient samples e.g. using Boule control material.

### Notes

- Prevent contamination of the reagents with dust, etc by use of the protective cap (delivered with the instrument).
- The cardboard box can be recycled after removal of the LPDE cubitainers. The LPDE container can be recycled after removal of any remaining reagent and a rinse with water.

### Ordering information and service

For international customers contact your local distributor for orders and support.

REF	Description	Packaging
1504122	Medonic M-series Diluent	20 L
1504123	Medonic M-series Lyse	5 L
1504128	Medonic M-series DualPack	1 L Lyse / 5 L Diluent
1504460	Medonic M-series Diluent, RFID	20 L
1504461	Medonic M-series Lyse, RFID	5 L
1504465	Medonic M-serie DualPack, RFID	1 L Lyse / 5 L Diluent

