

# Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22512-1K	2026-05-04
Low	22512-11	2026-05-04
Normal	22512-12	2026-05-04
High	22512-13	2026-05-04

Open vial stability 14 days <sup>1</sup>

MEDONIC		CONTROL	L
Medonic M-series/M32		LOT	22512-11
		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.28	0.12
MCV	fL	69.1	5.0
RDW	%	14.6	4.0
Hct	%	15.8	2.5
Plt	10 <sup>9</sup> /L	79	15
MPV	fL	11.7	1.5
WBC	10 <sup>9</sup> /L	3.2	0.3
Hgb	g/dL	5.8	0.3
MCH	pg	25.4	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 <sup>9</sup> /L	1.4	0.4
Mid#	10 <sup>9</sup> /L	0.2	0.3
Gran#	10 <sup>9</sup> /L	1.6	0.8
Lymph	%	43.6	7.0
Mid	%	5.7	5.0
Gran	%	50.7	7.0

CONTROL	N	
LOT	22512-12	
	Mean <sup>2</sup>	± Range <sup>3</sup>
	4.17	0.18
	82.6	5.0
	13.4	4.0
	34.4	3.0
	237	30
	11.2	1.5
	8.1	0.6
	12.0	0.4
	28.8	2.5
	34.9	3.0
	3.7	0.6
	0.5	0.5
	3.9	1.0
	45.5	5.0
	6.2	5.0
	48.3	5.0

CONTROL	H	
LOT	22512-13	
	Mean <sup>2</sup>	± Range <sup>3</sup>
	5.11	0.22
	92.0	5.0
	13.0	4.0
	47.0	3.5
	535	60
	11.5	1.5
	19.8	1.8
	16.1	0.5
	31.5	1.5
	34.3	3.0
	12.2	1.5
	1.1	1.0
	6.5	2.5
	61.3	5.0
	5.8	5.0
	32.9	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
  - hr Stabilnost otvorene bočice 14 dani (hr)
  - cs Stabilita otevřené ampulky 14 dní (cs)
  - da Åben rør stabilitet 14 dage (da)
  - et Stabiilsus avatud vialis 14 päevade (et)
  - fr Stabilité en flacon ouvert 14 jour (fr)
  - de Stabilität geöffneter Flaschen 14 tage (de)
  - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
  - it Stabilità della fiala aperta 14 giorni (it)
  - lv Atvērt flakonu stabilitāti 14 dienas (lv)
  - lt Stabilumas atidarius buteliuką 14 dienas (lt)
  - no Åpen rør/glass stabilitet 14 dager (no)
  - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
  - pt Estabilidade após abertura do frasco 14 dias (pt)
  - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
  - ru Стабильность вскрытого флакона 14 дней (ru)
  - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
  - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
  - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
  - es Estabilidad de la cápsula abierta 14 días (es)
  - sv Hållbarhet för öppnad flaska 14 dagar (sv)
  - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Sredno
  - Sredstvo
  - Průměr
  - Middelværdi
  - Keskmine
  - Moyenne
  - Mittelwert
  - Μέση
  - Media
  - Nozīmē
  - Vidurkis
  - Middelverdi
  - Średnia
  - Média
  - Medie
  - Среднее
  - Srednja vrednost
  - Priemer
  - Povprečno
  - Media
  - Medelvärde
  - Ortalama

- 3
- Обхват
  - Raspon
  - Rozsah
  - Område
  - Vahemik
  - Intervalle
  - Bereich
  - Αναμενόμενο εύρος α
  - Intervallo
  - Diapazons
  - Intervalas
  - Område
  - Zakres
  - Intervalo
  - Interval
  - диапазон
  - Opseg
  - Rozsah
  - Območje
  - Intervallo
  - Intervall
  - Aralığı

# For Medonic M-series / M32 with barcode reader

**M-Series:**

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

**M32:**

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series  
 Lot: 22512-11  
 Prod.: Con-Diff Low  
 ID: 2251211+  
 Exp.: 4-May-2026

Instr.: M-series  
 Lot: 22512-12  
 Prod.: Con-Diff Norm  
 ID: 2251212+  
 Exp.: 4-May-2026

Instr.: M-series  
 Lot: 22512-13  
 Prod.: Con-Diff High  
 ID: 2251213+  
 Exp.: 4-May-2026



3.

# Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22512-1K	2026-05-04
Low	22512-11	2026-05-04
Normal	22512-12	2026-05-04
High	22512-13	2026-05-04

Open vial stability 14 days <sup>1</sup>

## SWELAB

		CONTROL	L
		LOT	22512-11
Swelab Alfa/Alfa Plus		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.30	0.12
MCV	fL	70.4	5.0
RDW	%	13.9	4.0
Hct	%	16.2	2.5
Plt	10 <sup>9</sup> /L	81	15
MPV	fL	12.5	1.5
WBC	10 <sup>9</sup> /L	3.3	0.3
Hgb	g/dL	5.8	0.3
MCH	pg	25.2	1.5
MCHC	g/dL	35.8	3.0
Lymph#	10 <sup>9</sup> /L	1.6	0.4
Mid#	10 <sup>9</sup> /L	0.1	0.3
Gran#	10 <sup>9</sup> /L	1.6	0.8
Lymph	%	47.7	7.0
Mid	%	4.0	5.0
Gran	%	48.3	7.0

		CONTROL	N
		LOT	22512-12
		Mean <sup>2</sup>	± Range <sup>3</sup>
		4.17	0.18
		83.9	5.0
		12.3	4.0
		35.0	3.0
		243	30
		11.9	1.5
		8.2	0.6
		12.1	0.4
		29.0	1.5
		34.6	3.0
		4.1	0.6
		0.3	0.5
		3.8	1.0
		49.4	5.0
		4.2	5.0
		46.4	5.0

		CONTROL	H
		LOT	22512-13
		Mean <sup>2</sup>	± Range <sup>3</sup>
		5.11	0.22
		93.6	5.0
		12.1	4.0
		47.8	3.5
		550	60
		12.3	1.5
		19.8	1.8
		16.1	0.5
		31.5	1.5
		33.7	3.0
		12.8	1.5
		0.7	1.0
		6.3	2.5
		64.5	5.0
		3.5	5.0
		32.0	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
  - hr Stabilnost otvorene bočice 14 dani (hr)
  - cs Stabilita otevřené ampulky 14 dní (cs)
  - da Åben rør stabilitet 14 dage (da)
  - et Stabiilsus avatud vialis 14 päevade (et)
  - fr Stabilité en flacon ouvert 14 jour (fr)
  - de Stabilität geöffneter Flaschen 14 tage (de)
  - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
  - it Stabilità della fiala aperta 14 giorni (it)
  - lv Atvērt flakonu stabilitāti 14 dienas (lv)
  - lt Stabilumas atidarius buteliuką 14 dienas (lt)
  - no Åpen rør/glass stabilitet 14 dager (no)
  - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
  - pt Estabilidade após abertura do frasco 14 dias (pt)
  - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
  - ru Стабильность вскрытого флакона 14 дней (ru)
  - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
  - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
  - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
  - es Estabilidad de la cápsula abierta 14 días (es)
  - sv Hållbarhet för öppnad flaska 14 dagar (sv)
  - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
  - Sredstvo
  - Průměr
  - Middelverdi
  - Keskmine
  - Moyenne
  - Mittelwert
  - Μέση
  - Media
  - Nozīmē
  - Vidurkis
  - Middelverdi
  - Średnia
  - Média
  - Medie
  - Среднее
  - Srednja vrednost
  - Priemer
  - Povprečno
  - Media
  - Medelvärde
  - Ortalama

- 3
- Обхват
  - Raspon
  - Rozsah
  - Område
  - Vahemik
  - Intervalle
  - Bereich
  - Αναμενόμενο εύρος α
  - Intervallo
  - Diapazons
  - Intervalas
  - Område
  - Zakres
  - Intervalo
  - Interval
  - диапазон
  - Orseg
  - Rozsah
  - Območje
  - Intervalo
  - Intervall
  - Aralığı

## For Swelab Alfa / Alfa Plus with barcode reader

### Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

### Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa  
 Lot: 22512-11  
 Prod.: Con-Diff Low  
 ID: 2251211+  
 Exp.: 4-May-2026

Instr.: Swelab Alfa  
 Lot: 22512-12  
 Prod.: Con-Diff Norm  
 ID: 2251212+  
 Exp.: 4-May-2026

Instr.: Swelab Alfa  
 Lot: 22512-13  
 Prod.: Con-Diff High  
 ID: 2251213+  
 Exp.: 4-May-2026

