

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22603-1K	2026-08-03
Low	22603-11	2026-08-03
Normal	22603-12	2026-08-03
High	22603-13	2026-08-05

Open vial stability 14 days ¹

MEDONIC		CONTROL	L
Medonic M-series/M32		LOT	22603-11
		Mean ²	± Range ³
RBC	10 ¹² /L	2.29	0.12
MCV	fL	69.0	5.0
RDW	%	16.8	4.0
Hct	%	15.8	2.5
Plt	10 ⁹ /L	84	15
MPV	fL	11.3	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.8	0.3
MCH	pg	25.3	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	44.0	7.0
Mid	%	5.0	5.0
Gran	%	51.0	7.0

CONTROL	N	
LOT	22603-12	
	Mean ²	± Range ³
	4.23	0.18
	82.5	5.0
	13.5	4.0
	34.9	3.0
	239	30
	11.4	1.5
	8.5	0.6
	12.2	0.4
	28.8	2.5
	35.0	3.0
	4.0	0.6
	0.5	0.5
	4.0	1.0
	46.2	5.0
	6.2	5.0
	47.6	5.0

CONTROL	H	
LOT	22603-13	
	Mean ²	± Range ³
	5.15	0.22
	91.9	5.0
	13.0	4.0
	47.3	3.5
	537	60
	11.7	1.5
	20.8	1.8
	16.3	0.5
	31.7	1.5
	34.5	3.0
	12.8	1.5
	1.2	1.0
	6.8	2.5
	61.4	5.0
	5.7	5.0
	32.9	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dana (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelverdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervall
 - Aralıqı



www.boule.com

 **Boule**

 Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.:	M-series	Instr.:	M-series	Instr.:	M-series
Lot:	22603-11	Lot:	22603-12	Lot:	22603-13
Prod.:	Con-Diff Low	Prod.:	Con-Diff Norm	Prod.:	Con-Diff High
ID:	2260311+	ID:	2260312+	ID:	2260313+
Exp.:	3-Aug-2026	Exp.:	3-Aug-2026	Exp.:	5-Aug-2026

1:	 Bou l e C B A V 2	 Bou l e C B A V 2	 Bou l e C B A V 2
2:	 2260311+	 2260312+	 2260313+
3:	 26080300001153030000	 26080300002153030000	 26080500003153030000
4:	 10004029800000000000	 10004029700000000000	 10004029300000000000
5:	 22912690500580315825	 42318825501220434930	 51522919501630547335
6:	 08415113152531536730	 23930114152881535030	 53760117153171534530
7:	 03403015040020301708	 08506040060050504010	 20818128150121006825
8:	 16840440700505051070	 13540462500625047650	 13040614500575032950
9:	 11119422521118226312	 11114900993234169806	 11117919122370543872

3.

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22603-1K	2026-08-03
Low	22603-11	2026-08-03
Normal	22603-12	2026-08-03
High	22603-13	2026-08-05

Open vial stability 14 days ¹

SWELAB

CONTROL		L	
LOT		22603-11	
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.30	0.12
MCV	fL	70.1	5.0
RDW	%	15.8	4.0
Hct	%	16.1	2.5
Plt	10 ⁹ /L	83	15
MPV	fL	11.8	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.9	0.3
MCH	pg	25.7	1.5
MCHC	g/dL	36.6	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	46.9	7.0
Mid	%	3.3	5.0
Gran	%	49.8	7.0

CONTROL		N	
LOT		22603-12	
		Mean ²	± Range ³
		4.23	0.18
		83.8	5.0
		12.5	4.0
		35.4	3.0
		242	30
		12.1	1.5
		8.5	0.6
		12.2	0.4
		28.8	1.5
		34.5	3.0
		4.2	0.6
		0.3	0.5
		4.0	1.0
		49.6	5.0
		3.5	5.0
		46.9	5.0

CONTROL		H	
LOT		22603-13	
		Mean ²	± Range ³
		5.12	0.22
		93.6	5.0
		12.1	4.0
		47.9	3.5
		546	60
		12.4	1.5
		20.8	1.8
		16.2	0.5
		31.6	1.5
		33.8	3.0
		13.5	1.5
		0.7	1.0
		6.6	2.5
		64.5	5.0
		3.6	5.0
		31.9	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilitnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelværdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Orseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı



www.boule.com

 **Boule**



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 22603-11
Prod.: Con-Diff Low
ID: 2260311+
Exp.: 3-Aug-2026

Instr.: Swelab Alfa
Lot: 22603-12
Prod.: Con-Diff Norm
ID: 2260312+
Exp.: 3-Aug-2026

Instr.: Swelab Alfa
Lot: 22603-13
Prod.: Con-Diff High
ID: 2260313+
Exp.: 5-Aug-2026

