

Certificate of Analysis


Product	MEDONIC MEDILYSE II-Vet (5 L)
Product No	1503975
Lot No	1801-878
Manufacturing Date	2018-01-15
Expiration Date	2021-01-15


Test	Specification
pH	7.10 ± 0.15
Osmolality (mOsm/kg)	210 ± 15
Conductivity (mS/cm)	12.45 ± 0.35
CBC Performance *	Passes test
Particles (PLT/RBC channel)	≤ 5
Total bacterial growth (CFU/100ml)	≤ 5

* Performance test using calibrated Blood cell counter at Boule Medical reference laboratory. All parameters except mid cells are tested. Accepted reference 95 % to 105 % of expected control blood sample values.

Actual values meet specifications.

Date: 2018-01-31


Jennifer Abrahamsson
Quality Assurance



This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

