

Certificate of Analysis

| | |
|---------------------------|-----------------------------|
| Product | MEDONIC M51-D DILUENT (20L) |
| Product No | 1504510 |
| Lot No | 20181231 |
| Manufacturing Date | 2018-12-31 |
| Expiration Date | 2020-12-30 |


| Test | Specification |
|---|----------------------|
| pH | 7,4 – 7,8 |
| Osmolality (mOsm/kg) | 300 – 320 |
| Conductivity (mS/cm) | 15,5 – 16,5 |
| Particles PLT/RBC channel ($\times 10^9/L$) | ≤ 5 |

Actual values meet specification.

Date: 2019-03-12



Lubna Elabbas
Quality Assurance



This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

