

CERTIFICATE OF ANALYSIS

Product: Swelab™AlfaDiluent

Lot No: 2011-002

Product Part #: 501-208


Exp. Date: 2023-10-31

Manufacturing Date: 2020-11-13

Quality Testing Date: 2020-11-20

Parameter	Specification	Result	Pass/Fail
pH (units):	6.50 – 6.90	6.78	Pass
Osmolality (mOsm/kg):	294 – 312	298	Pass
Conductivity (mS/cm):	16.98 – 17.93	17.33	Pass
Background (x10 ³ /	0.0 – 8.0	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe
Quality Assurance Manager

12/4/2020

Date

