

CERTIFICATE OF ANALYSIS

Product: Medonic™ M-Series Lyse

Lot No: 2106-002

Product Part #: 501-212

Exp. Date: 2024-05-31

Manufacturing Date: 2021-06-04

Quality Testing Date: 2021-06-08

Parameter	Specification	Result	Pass/Fail
pH (units):	6.80 – 7.20	7.13	Pass
Osmolality (mOsm/kg):	220 – 240	233	Pass
Conductivity (mS/cm):	14.00 – 15.00	14.20	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Louise Esteva
Regulatory and Compliance Specialist

06-23-2021

Date

