

# Certificate of Analysis


|                           |                       |
|---------------------------|-----------------------|
| <b>Product</b>            | SWELAB ALFALYSE (5 L) |
| <b>Product No</b>         | 1504125               |
| <b>Lot No</b>             | 2109-568              |
| <b>Manufacturing Date</b> | 2021-09-01            |
| <b>Expiration Date</b>    | 2024-09-01            |

| Test                    | Specification |
|-------------------------|---------------|
| pH                      | 7,80 – 8,30   |
| Osmolality (mOsm/kg)    | 295 – 330     |
| Conductivity (mS/cm)    | 17,75 – 18,25 |
| CBC Performance *       | Passes test   |
| Particles (WBC channel) | ≤ 0,2         |


\* Performance test using calibrated Blood cell counter at Boule Medical reference laboratory. All parameters except mid cells are tested. Accepted reference 95 % to 105 % of expected control blood sample values.

Actual value meets specification.

Date: 2021-09-09



Lubna Elabbas  
Quality Assurance



This product is intended for Laboratory use only unless otherwise specifically stated.  
All information in this document is property of Boule Medical AB.

