

CERTIFICATE OF ANALYSIS

Product: Medonic™ M-Series Diluent

Lot No: 2107-003

Product Part #: 501-211

Exp. Date: 2024-06-30

Manufacturing Date: 2021-07-29

Quality Testing Date: 2021-08-02

| Parameter | Specification | Result | Pass/Fail |
|---|---------------|--------|-----------|
| pH (units): | 6.90 – 7.20 | 7.13 | Pass |
| Osmolality (mOsm/kg): | 290 – 310 | 292 | Pass |
| Conductivity (mS/cm): | 17.50 – 18.00 | 17.58 | Pass |
| Background Test (x10 ³ /μL): | 0.0 – 5.0 | 0.5 | Pass |

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe
Quality and Regulatory Manager

8/4/2021

Date

