

CERTIFICATE OF ANALYSIS

Product: Swelab™ AlfaDiluent

Lot No: 2109-001

Product Part #: 501-208


Exp. Date: 2024-08-31

Manufacturing Date: 2021-09-28

Quality Testing Date: 2021-09-29

Parameter	Specification	Result	Pass/Fail
pH (units):	6.50 – 6.90	6.74	Pass
Osmolality (mOsm/kg):	294 – 312	299	Pass
Conductivity (mS/cm):	16.98 – 17.93	17.35	Pass
Background (x10 ³ /μL):	0.0 – 8.0	0.5	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe
Quality Assurance Manager

10/04/2024

Date

