

CERTIFICATE OF ANALYSIS

Product: Swelab™ Alfalyse, 5L**Lot No:** 2107-001**Product Part #:** 501-207**Exp. Date:** 2024-06-30**Manufacturing Date:** 2021-07-14**Quality Testing Date:** 2021-07-16

Parameter	Specification	Result	Pass/Fail
pH (units):	7.54 – 8.30	8.27	Pass
Osmolality (mOsm/kg):	295 – 330	318	Pass
Conductivity (mS/cm):	18.00– 18.50	18.34	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe
Quality Assurance Manager



Date

