

# Certificate of Analysis

|                           |                                   |
|---------------------------|-----------------------------------|
| <b>Product</b>            | MEDONIC M-SERIES LYSE, RFID (5 L) |
| <b>Product No</b>         | 1504461                           |
| <b>Lot No</b>             | 2110-598                          |
| <b>Manufacturing Date</b> | 2021-10-20                        |
| <b>Expiration Date</b>    | 2024-10-20                        |

| Test                    | Specification |
|-------------------------|---------------|
| pH                      | 6,80 – 7,20   |
| Osmolality (mOsm/kg)    | 220 – 240     |
| Conductivity (mS/cm)    | 14,00 – 15,00 |
| CBC Performance test*   | Passes test   |
| Particles (WBC channel) | $\leq 0,2$    |


\* Performance test using calibrated Blood cell counter at Boule Medical reference laboratory. All parameters except mid cells are tested. Accepted reference 95 % to 105 % of expected control blood sample values.

Actual value meets specification.

Date: 2021-10-27



Lubna Elabbas  
Quality Assurance



This product is intended for Laboratory use only unless otherwise specifically stated.  
All information in this document is property of Boule Medical AB.

