

## CERTIFICATE OF ANALYSIS

**Product:** Medonic™ M-Series Lyse

**Lot No:** 2112-002

**Product Part #:** 501-212


**Exp. Date:** 2024-11-30

**Manufacturing Date:** 12/14/2021

**Quality Testing Date:** 12/23/2021

Parameter	Specification	Result	Pass/Fail
pH (units):	6.80 – 7.20	7.10	Pass
Osmolality (mOsm/kg):	220 – 240	229	Pass
Conductivity (mS/cm):	14.00 – 15.00	14.14	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

  
\_\_\_\_\_  
Jeanne Roscoe  
Regulatory and Compliance Manager

Date

  
\_\_\_\_\_  
1/21/22