

CERTIFICATE OF ANALYSIS

Product: Swelab™ AlfaDiluent RFID

Lot No: 2109-001

Product Part #: 501-465


Exp. Date: 2024-08-31

Manufacturing Date: 2021-09-28

Quality Testing Date: 2021-09-29

| Parameter | Specification | Result | Pass/Fail |
|------------------------------------|---------------|--------|-----------|
| pH (units): | 6.50 – 6.90 | 6.74 | Pass |
| Osmolality (mOsm/kg): | 294 – 312 | 297 | Pass |
| Conductivity (mS/cm): | 16.98 – 17.93 | 17.35 | Pass |
| Background (x10 ³ /μL): | 0.0 – 8.0 | 0.5 | Pass |

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe
Quality Assurance Manager

10/04/2021

Date

