

CERTIFICATE OF ANALYSIS

Product: Swelab™ AlfaLyse RFID

Lot No: 2203-002

Product Part #: 501-466


Exp. Date: 2025-02-28

Manufacturing Date: 2022-03-08

Quality Testing Date: 2022-03-11

Parameter	Specification	Result	Pass/Fail
pH (units):	7.54 – 8.30	8.28	Pass
Osmolality (mOsm/kg):	295 – 330	316	Pass
Conductivity (mS/cm):	18.00 – 18.50	18.32	Pass
Functional Test	Pass/Fail	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe
Quality Assurance Manager

4/20/22

Date

