

CERTIFICATE OF ANALYSIS

Product: Swelab™ AlfaLyse

Lot No: 1669

Product Part #: 501-210

Exp. Date: 2025-04-30

Manufacturing Date: 05/18/2022

Quality Testing Date: 05/23/2022

Parameter	Specification	Result	Pass/Fail
pH (units):	7.54 – 8.30	8.27	Pass
Osmolality (mOsm/kg):	295 – 330	319	Pass
Conductivity (mS/cm):	18.00 – 18.50	18.28	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe
Quality Assurance Manager

5/25/2022

Date

