


# Certificate of Analysis


|                           |                                 |
|---------------------------|---------------------------------|
| <b>Product</b>            | MEDONIC M-SERIES DILUENT (20 L) |
| <b>Product No</b>         | 1504122                         |
| <b>Lot No</b>             | 2205-715                        |
| <b>Manufacturing Date</b> | 2022-05-17                      |
| <b>Expiration Date</b>    | 2025-05-17                      |

| <b>Test</b>                 | <b>Specification</b> |
|-----------------------------|----------------------|
| pH                          | 6,90 – 7,10          |
| Osmolality (mOsm/kg)        | 290 – 310            |
| Conductivity (mS/cm)        | 17,50 – 18,00        |
| Particles (PLT/RBC channel) | ≤ 5                  |

Actual value meets specification.



Date: 2022-05-24

  
*Ulrika Miefalk*  
Ulrika Miefalk  
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.  
All information in this document is property of Boule Medical AB.

