

## CERTIFICATE OF ANALYSIS

**Product:** Medonic™ M-Series Lyse

**Lot No:** 2206-004

**Product Part #:** 501-212

**Exp. Date:** 2025-05-31

**Manufacturing Date:** 2022-06-10


**Quality Testing Date:** 2022-06-15

| Parameter             | Specification | Result | Pass/Fail |
|-----------------------|---------------|--------|-----------|
| pH (units):           | 6.80 – 7.20   | 7.12   | Pass      |
| Osmolality (mOsm/kg): | 220 – 240     | 233    | Pass      |
| Conductivity (mS/cm): | 14.00 – 15.00 | 14.18  | Pass      |
| Functional Test       | Pass          | Pass   | Pass      |

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe  
Regulatory and Compliance Manager



Date

