

Boule Con-5Diff A1

CONTROL

IVD

LOT



Tri-Level	72207-1K	2022-11-05
Low	72207-11	2022-11-05
Normal	72207-12	2022-11-05
High	72207-13	2022-11-05

www.medonic.se
www.swelab.se
Support/downloads

Open vial stability: 14 days ¹

Medonic M51 / Swelab Lumi

Parameter ²	Unit ³	Low		Normal		High	
		Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵
WBC	10 ⁹ /L	3.52	0.50	8.37	1.00	19.34	2.50
Neu%	%	50.5	16.0	55.5	8.0	62.8	8.0
Lym%	%	37.9	21.0	29.8	8.0	19.7	12.0
Mon%	%	5.2	5.0	5.7	5.0	6.1	6.0
Eos%	%	6.4	5.0	9.0	6.0	11.4	8.0
Bas%	%	4.6	10.0	3.6	10.0	2.7	10.0
Neu#	10 ⁹ /L	1.78	0.50	4.66	0.70	12.15	2.00
Lym#	10 ⁹ /L	1.33	0.70	2.49	0.80	3.81	2.00
Mon#	10 ⁹ /L	0.18	0.30	0.47	0.50	1.18	1.10
Eos#	10 ⁹ /L	0.22	0.20	0.75	0.60	2.20	1.40
Bas#	10 ⁹ /L	0.16	0.40	0.30	0.80	0.52	1.80
RBC	10 ¹² /L	2.36	0.24	4.62	0.30	5.52	0.50
HGB	g/L	60	6	133	6	173	8
HCT	%	19.4	2.0	41.8	3.0	54.2	4.0
MCV	fL	82.0	5.0	90.5	5.0	98.2	6.0
MCH	pg	25.4	2.5	28.7	2.5	31.4	2.5
MCHC	g/L	309	30	318	30	319	30
RDW-CV	%	15.2	3.0	13.4	3.0	13.0	3.0
PLT	10 ⁹ /L	56	20	254	40	455	60
MPV	fL	8.7	3.0	8.6	3.0	8.8	3.0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.

<p>1 bg Трайност на отворен флакон: 14 дни cs Stabilita otevřené ampulky: 14 dny da Åben rør stabilitet: 14 dage de Stabilität geöffneter Flaschen: 14 tage el Διατήρηση δείγματος μετά την αποσφράγιση: 14 ημέρα es Estabilidad de la cápsula abierta: 14 días et Stabiilsus avatud vialis: 14 päevade arv fr Stabilité en flacon ouvert: 14 jour hr Stabilnost otvorene bočice: 14 dani it Stabilità della fiala aperta: 14 giorni lt Stabilumas atidarius buteliuka: 14 dienų lv Atvērt flakonu stabilitāte: 14 dienas no Åpen rør/glass stabilitet: 14 dager pl Trwałość po otwarciu fioleki: 14 Liczba dni pt Estabilidade após abertura do frasco: 14 dias ro Termenul de valabilitate al fiolele desfăcute: 14 zile ru Стабильность вскрытого флакона: 14 дней sk Stabilita otvorenej injekčnej liekovky: 14 dní sl Stabilnost odprte stekleničke: 14 dni sr Proizvod stabilan nakon otvaranja pakovanja: 14 dana sv Hållbarhet för öppnad flaska: 14 dagar tr Acik sise stabilitesi aün: 14 gün</p>	<p>2 Параметър Parametr Parameter Parameter Παράμετρος Parámetro Parameeter Paramètre Parametar Parametro Parametras Parametrs Parameter Kryterium Parâmetro Parametru Параметр Parameter Parameter Parameter Parameter Parametar Parameter Parametre</p>	<p>3 Единица Jednotka Enhed Einheit Einheit Μονάδα Unidad Ühik Unité Jedinica Parametar Unità Vienetas Vienība Enhet Jednostka Unidade Unitate Единица измерения Jednotka Enota Jedinica Enhet Birim</p>	<p>4 Средно Průměr Middelværdi Mittelwert Μέση Media Keskmine Moyenne Sredstvo Media Vidurkis Nozīmē Middelverdi Średnia Média Medie Среднее Priemer Povprečno Srednja vrednost Medelvärde Ortalama</p>	<p>5 Обхват Rozsah Område Bereich Αναμενόμενο εύρος α Intervalo Vahemik Intervalle Raspon Intervalo Intervalas Diapazons Område Zakres Intervalo Interval Rozsah Območje Opseg Intervall Aralıđı</p>
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

31340-8

1504546