

CERTIFICATE OF ANALYSIS

Product: Swelab™AlfaDiluent

Lot No: 2209-004

Product Part #: 501-208

Exp. Date: 2025-08-31

Manufacturing Date: 2022-09-29

Quality Testing Date: 2022-10-03

| Parameter | Specification | Result | Pass/Fail |
|------------------------------------|---------------|--------|-----------|
| pH (units): | 6.50 – 6.90 | 6.84 | Pass |
| Osmolality (mOsm/kg): | 294 – 312 | 297 | Pass |
| Conductivity (mS/cm): | 16.98 – 17.93 | 17.18 | Pass |
| Background (x10 ³ /μL): | 0.0 – 8.0 | 0.0 | Pass |

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Rescoe
Quality Assurance Manager

10/12/22

Date

