

CERTIFICATE OF ANALYSIS

Product: Swelab™ AlfaLyse RFID

Lot No: 2109-002

Product Part #: 501-466


Exp. Date: 2024-08-31

Manufacturing Date: 2021-09-29

Quality Testing Date: 2021-10-04

Parameter	Specification	Result	Pass/Fail
pH (units):	7.54 – 8.30	8.26	Pass
Osmolality (mOsm/kg):	295 – 330	316	Pass
Conductivity (mS/cm):	18.00 – 18.50	18.35	Pass
Functional Test	Pass/Fail	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe
Quality Assurance Manager

10/5/2021

Date

