

Certificate of Analysis

Product	SWELAB LUMI-L2 LYSE (500 mL)
Product No	1504516
Lot No	2022090401
Manufacturing Date	2022-09-04
Expiration Date	2024-09-03

Test	Specification
pH	4,0 – 4,4
CBC Performance test*	Passes test
Background WBC ($\times 10^9/L$)	$\leq 0,2$

*Performance testing at a reference laboratory. WBC, HGB and WBC scattergrams tested and approved.

Actual values meet specification.

Date: 2022-10-18



Sofia Karlström
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

