

Certificate of Analysis


| | |
|---------------------------|-----------------------------------|
| Product | MEDONIC M-SERIES LYSE, RFID (5 L) |
| Product No | 1504461 |
| Lot No | 2301-851 |
| Manufacturing Date | 2023-01-12 |
| Expiration Date | 2026-01-12 |

| Test | Specification |
|-------------------------|----------------------|
| pH | 6,80 – 7,20 |
| Osmolality (mOsm/kg) | 220 – 240 |
| Conductivity (mS/cm) | 14,00 – 15,00 |
| CBC Performance test* | Passes test |
| Particles (WBC channel) | ≤ 0,2 |

* Performance test using calibrated Blood cell counter at Boule Medical reference laboratory. All parameters except mid cells are tested. Accepted reference 95 % to 105 % of expected control blood sample values.

Actual value meets specification.

Date: 2023-01-18



Ulrika Miefalk
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

