

CERTIFICATE OF ANALYSIS

Product: Swelab™Alfalyse, 5L

Lot No: 2205-003

Product Part #: 501-207

Exp. Date: 2025-04-30

Manufacturing Date: 2022-05-31


Quality Testing Date: 2022-06-03

Parameter	Specification	Result	Pass/Fail
pH (units):	7.54 – 8.30	8.29	Pass
Osmolality (mOsm/kg):	295 – 330	318	Pass
Conductivity (mS/cm):	18.00– 18.50	18.32	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Jeanne Roscoe
Quality Assurance Manager



Date

