

## CERTIFICATE OF ANALYSIS

**Product:** Swelab™AlfaDiluent

**Lot No:** 2301-002

**Product Part #:** 501-208

**Exp. Date:** 2025-12-31

**Manufacturing Date:** 01-24-2023

**Quality Testing Date:** 01-26-2023

Parameter	Specification	Result	Pass/Fail
pH (units):	6.50 – 6.90	6.82	Pass
Osmolality (mOsm/kg):	294 – 312	299	Pass
Conductivity (mS/cm):	16.98 – 17.93	17.30	Pass
Background (x10 <sup>3</sup> /μL):	0.0 – 8.0	0.5	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

  
\_\_\_\_\_  
Quality Assurance

2/3/2023  
Date

