

CERTIFICATE OF ANALYSIS

Product: Quintus 5-Part Lyse**Lot No:** 2301-301**Product Part #:** 501-291**Exp. Date:** 2024-07-31**Manufacturing Date:** 2023-01-09**Quality Testing Date:** 2023-01-11

Parameter	Specification	Result	Pass/Fail
pH (units):	2.7 – 3.0	3.0	Pass
Osmolality (mOsm/kg):	35 – 44	37	Pass
Conductivity (µS/cm):	1,117 – 1,541	1,253	Pass


The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Quality & Regulatory Department

2/16/23

Date



CERTIFICATE OF ANALYSIS

Product: Quintus 5-Part Diluent

Lot No: 2301-201

Product Part #: 501-292

Exp. Date: 2024-12-31

Manufacturing Date: 2023-01-20

Quality Testing Date: 2023-02-10

Parameter	Specification	Result	Pass/Fail
pH (units):	6.7 – 7.1	7.1	Pass
Osmolality (mOsm/kg):	320 – 350	325	Pass
Conductivity (mS/cm):	18.50 – 19.50	18.68	Pass
Background ($\times 10^3/\mu\text{L}$)	0 – 5	3	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.


Quality & Regulatory Department

Date 2/16/23

