

Certificate of Analysis

Product	SWELAB LUMI-L2 LYSE (500 mL)
Product No	1504516
Lot No	2022111901
Manufacturing Date	2022-11-19
Expiration Date	2024-11-18

Test	Specification
pH	4,0 – 4,4
CBC Performance test*	Passes test
Background WBC (x 10 ⁹ /L)	≤ 0,2

*Performance testing at a reference laboratory. WBC, HGB and WBC scattergrams tested and approved.

Actual values meet specification.

Date: 2023-03-02



Ulrika Miefalk
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

