

CERTIFICATE OF ANALYSIS

Product: Medonic™ M-Series Diluent

Lot No: 2303-003

Product Part #: 501-211

Exp. Date: 2026-02-28

Manufacturing Date: 2023-03-28

Quality Testing Date: 2023-03-31

| Parameter | Specification | Result | Pass/Fail |
|---|---------------|--------|-----------|
| pH (units): | 6.90 – 7.20 | 7.09 | Pass |
| Osmolality (mOsm/kg): | 290 – 310 | 291 | Pass |
| Conductivity (mS/cm): | 17.50 – 18.00 | 17.55 | Pass |
| Background Test (x10 ³ /μL): | 0.0 – 5.0 | 0.5 | Pass |

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Quality and Regulatory

4/10/23

Date

