


Certificate of Analysis

| | |
|---------------------------|---|
| Product | Boule Micropipette Plastic EDTA 10x100 pcs Boule Micropipette Plastic EDTA 1x100 pcs |
| Product No | 1070030 1070039 |
| Lot No | 230145D5 |
| Manufacturing Date | 2023-03-09 |
| Expiration Date | 2026-02-28 |

| Parameter | Specification |
|-------------------|------------------------------|
| Volume 20 μ L | $CV \leq 1,0\%$ |
| WBC background | $\leq 0,2 \times 10^9/L$ |
| RBC background | $\leq 0,03 \times 10^{12}/L$ |
| HGB background | $\leq 0,2 \text{ x g/dL}$ |
| PLT background | $\leq 10 \times 10^9/L$ |

Actual value meets specification.

Date: 2023-03-20



Caroline Sparf
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

