

## CERTIFICATE OF ANALYSIS

**Product:** Swelab™ AlfaDiluent

**Lot No:** 2304-003

**Product Part #:** 501-208

**Exp. Date:** 2026-03-31

**Manufacturing Date:** 04-28-2023

**Quality Testing Date:** 05-02-2023

Parameter	Specification	Result	Pass/Fail
pH (units):	6.50 – 6.90	6.83	Pass
Osmolality (mOsm/kg):	294 – 312	301	Pass
Conductivity (mS/cm):	16.98 – 17.93	17.27	Pass
Background (x10 <sup>3</sup> /μL):	0.0 – 8.0	1.0	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

  
\_\_\_\_\_  
Quality Assurance

5/16/23  
\_\_\_\_\_  
Date