

Certificate of Analysis

Product H50V SHEATH, RFID
Product No 1440107 (5 L)
Lot No 23041412M2
Manufacturing Date 2023-04-26
Expiration Date 2025-04-25

Test	Specification
Appearance	Sheath should be transparent liquid without precipitation, particles or floccules
pH	8.55 - 8.65
Osmolality (mOsm/kg)	20 - 60
Blanks	WBC $\leq$ 0.2 $\times$ 10 ⁹ /L

Actual values meet specifications.

Date: 2023-06-29



Pernilla Hammarstedt Rydberg
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

