

## CERTIFICATE OF ANALYSIS

**Product:** Swelab™Alfalyse, 5L

**Lot No:** 1879

**Product Part #:** 501-207

**Exp. Date:** 2026-06-30

**Manufacturing Date:** 2023-07-18

**Quality Testing Date:** 2023-07-24

Parameter	Specification	Result	Pass/Fail
pH (units):	7.54 – 8.30	8.28	Pass
Osmolality (mOsm/kg):	295 – 330	320	Pass
Conductivity (mS/cm):	18.00– 18.50	18.26	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

  
Quality Assurance

07/26/23  
Date



