

# Certificate of Analysis

**Product** Quintus 5-part Diluent 20L  
**Product No** 1504291  
**Lot No** 2307-935  
**Manufacturing Date** 2023-07-03  
**Expiration Date** 2025-07-02

Test	Specification
pH	6,8 – 7,2
Osmolality (mOsm/kg)	315 – 350
Conductivity (mS/cm)	18,8 – 19,2
Background PLT ( x 10 <sup>9</sup> /L)	≤ 5

Actual values meet specification.

Date: 2023-08-17



*Sofia Karlström*  
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.  
All information in this document is property of Boule Medical AB.

