

CERTIFICATE OF ANALYSIS

Product: Medonic™ M-Series Lyse

Lot No: 2308-002

Product Part #: 501-212


Exp. Date: 2026-07-31

Manufacturing Date: 2023-08-21

Quality Testing Date: 2023-08-28

| Parameter | Specification | Result | Pass/Fail |
|-----------------------|---------------|--------|-----------|
| pH (units): | 6.80 – 7.20 | 7.14 | Pass |
| Osmolality (mOsm/kg): | 220 – 240 | 233 | Pass |
| Conductivity (mS/cm): | 14.00 – 15.00 | 14.16 | Pass |
| Functional Test | Pass | Pass | Pass |

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.


Regulatory and Compliance

08/29/23.
Date

