

# CERTIFICATE OF ANALYSIS

**Product:** Boule Cleaning Kit

**Kit Lot No:** 1899

**Product Part #:** 501-036


**Kit Exp. Date:** 2025-08-31

<b>Part Number:</b> 503-085			
<b>Boule Detergent Cleaner Lot #</b> 1900			
<b>Manufacturing Date:</b> 2023-08-08		<b>Quality Testing Date:</b> 2023-08-14	
Parameter	Specification	Result	Pass/Fail
pH (units):	7.60 – 8.40	8.29	Pass

<b>Part Number:</b> 503-086			
<b>Boule Hypochlorite (2%) Cleaner Lot #</b> 1901			
<b>Manufacturing Date:</b> 2023-08-08		<b>Quality Testing Date:</b> 2023-08-10	
Parameter	Specification	Result	Pass/Fail
pH (units):	11.80 – 13.00	12.98	Pass

<b>Part Number:</b> 503-087			
<b>Boule Enzymatic Cleaner Lot #</b> 1902			
<b>Manufacturing Date:</b> 2023-08-10		<b>Quality Testing Date:</b> 2023-08-18	
Parameter	Specification	Result	Pass/Fail
pH (units):	7.60 – 8.40	8.10	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

  
Quality and Regulatory

08/18/23.  
Date

