

CERTIFICATE OF ANALYSIS

Product: Quintus 5-Part Stopper

Lot No: 2310-201

Product Part #: 501-350


Exp. Date: 2026-09-30

Manufacturing Date: 2023-10-09

Quality Testing Date: 2023-10-11

Parameter	Specification	Result	Pass/Fail
pH (units):	10.10 – 11.20	11.08	Pass
Osmolality (mOsm/kg):	830 – 860	843	Pass
Conductivity (mS/cm):	42.7 – 43.3	43.1	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Quality Assurance

10/26/23

Date

