

# Certificate of Analysis

<b>Product</b>	EXIGO EOS REAGENT, RFID (1.9 L)
<b>Product No</b>	1504502
<b>Lot No</b>	2401-136
<b>Manufacturing Date</b>	2024-01-15
<b>Expiration Date</b>	2026-01-14

<b>Test</b>	<b>Specification</b>
pH	8.70 – 9.30
Osmolality (mOsm/kg)	125 - 145
Conductivity (mS/cm)	6.80 – 7.60
EOS Performance *	Passes test
Blank EOS	≤ 0.1

\* Performance test using calibrated Blood cell counter at Boule Medical reference laboratory.

Actual values meet specifications.

Date: 2024-01-17



Ulrika Miefalk  
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.  
All information in this document is property of Boule Medical AB.

