

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	---	---
Low	22401-41	2024-05-30
Normal	22401-42	2024-05-30
High	22401-43	2024-06-04

MEDONIC

CONTROL		L	
LOT		22401-41	
Medonic M-series/M32		Mean ²	± Range ³
RBC	10 ¹² /L	2.23	0.12
MCV	fL	66.1	5.0
RDW	%	14.2	4.0
Hct	%	14.7	2.5
Plt	10 ⁹ /L	79	15
MPV	fL	11.1	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	25.1	1.5
MCHC	g/dL	38.1	3.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.6	0.8
Lymph	%	45.8	7.0
Mid	%	5.0	5.0
Gran	%	49.2	7.0

CONTROL		N	
LOT		22401-42	
		Mean ²	± Range ³
		4.18	0.18
		78.5	5.0
		14.2	4.0
		32.8	3.0
		225	30
		11.4	1.5
		8.3	0.6
		11.7	0.4
		28.0	1.5
		35.7	3.0
		3.9	0.6
		0.4	0.5
		4.0	1.0
		46.7	5.0
		5.0	5.0
		48.3	5.0

CONTROL		H	
LOT		22401-43	
		Mean ²	± Range ³
		5.08	0.22
		91.7	5.0
		13.6	4.0
		46.6	3.5
		513	60
		11.9	1.5
		20.6	1.8
		16.2	0.5
		31.9	1.5
		34.8	3.0
		12.7	1.5
		1.0	1.0
		6.9	2.5
		61.5	5.0
		5.0	5.0
		33.5	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilnost otvorene bočice 14 dni (hr)
cs	Stabilita otevřených ampulky 14 dní (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud viaalis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διάρκεια δειγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiolek 14 Liczba dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvorene injekčneje liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 días (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

Средно
Sredstvo
Průměr
Middelværdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nozīmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3

Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Aralığı

For Medonic M-series / M32 with barcode reader

- M-Series:**
1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

- M32:**
1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9
- | | | | |
|---------|--------------|---------|---------------|
| Instr.: | M-series | Instr.: | M-series |
| Lot: | 22401-41 | Lot: | 22401-42 |
| Prod.: | Con-Diff Low | Prod.: | Con-Diff Norm |
| ID: | 2240141+ | ID: | 2240142+ |
| Exp.: | 30-May-2024 | Exp.: | 30-May-2024 |

- | | |
|---------|---------------|
| Instr.: | M-series |
| Lot: | 22401-43 |
| Prod.: | Con-Diff High |
| ID: | 2240143+ |
| Exp.: | 4-Jun-2024 |

1:	 Boul e CBAV2	 Boul e CBAV2	 Boul e CBAV2
2:	 2240141+	 2240142+	 2240143+
3:	 24053000001153030000	 24053000002153030000	 24060400003153030000
4:	 10004031200000000000	 10004031100000000000	 10004031100000000000
5:	 22312661500560314725	 41818785501170432830	 50822917501620546635
6:	 07915111152511538130	 22530114152801535730	 51360119153191534830
7:	 03303015040020301608	 08306039060040504010	 20618127150101006925
8:	 14240458700505049270	 14240467500505048350	 13640615500505033550
9:	 11116774039063832776	 11116834562408122259	 11118881311530054881

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	---	---
Low	22401-41	2024-05-30
Normal	22401-42	2024-05-30
High	22401-43	2024-06-04

SWELAB

CONTROL		L	
LOT		22401-41	
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.23	0.12
MCV	fL	68.4	5.0
RDW	%	14.0	4.0
Hct	%	15.3	2.5
Plt	10 ⁹ /L	81	15
MPV	fL	12.3	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	25.1	1.5
MCHC	g/dL	36.6	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	48.0	7.0
Mid	%	3.1	5.0
Gran	%	48.9	7.0

CONTROL		N	
LOT		22401-42	
		Mean ²	± Range ³
		4.22	0.18
		80.3	5.0
		13.2	4.0
		33.9	3.0
		233	30
		12.5	1.5
		8.4	0.6
		11.9	0.4
		28.2	1.5
		35.1	3.0
		4.2	0.6
		0.3	0.5
		3.9	1.0
		50.1	5.0
		3.3	5.0
		46.6	5.0

CONTROL		H	
LOT		22401-43	
		Mean ²	± Range ³
		5.10	0.22
		93.5	5.0
		12.5	4.0
		47.7	3.5
		534	60
		12.9	1.5
		20.6	1.8
		16.4	0.5
		32.2	1.5
		34.4	3.0
		13.3	1.5
		0.6	1.0
		6.7	2.5
		64.5	5.0
		3.0	5.0
		32.5	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilitet otvorene bočice 14 dni (hr)
cs	Stabilita otevřených ampulky 14 dní (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud vialis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διάρθρωση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiolek 14 dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvoreneje injekčneje liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 días (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

Sredstvo
Průměr
Middelverdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nožmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3

Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Aralığı

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 22401-41
Prod.: Con-Diff Low
ID: 2240141+
Exp.: 30-May-2024

Instr.: Swelab Alfa
Lot: 22401-42
Prod.: Con-Diff Norm
ID: 2240142+
Exp.: 30-May-2024

Instr.: Swelab Alfa
Lot: 22401-43
Prod.: Con-Diff High
ID: 2240143+
Exp.: 4-Jun-2024

