

## CERTIFICATE OF ANALYSIS

**Product:** Quintus 5-Part Diluent

**Lot No:** 2401-201

**Product Part #:** 501-292

**Exp. Date:** 2025-12-31

**Manufacturing Date:** 2024-01-29


**Quality Testing Date:** 2024-02-27

Parameter	Specification	Result	Pass/Fail
pH (units):	6.7 – 7.1	7.1	Pass
Osmolality (mOsm/kg):	320 – 350	328	Pass
Conductivity (mS/cm):	18.50 – 19.50	18.66	Pass
Background ( x10 <sup>3</sup> /μL)	0 – 5	4	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Quality and Regulatory



Date

