

Boule Con-5Diff A1

CONTROL

IVD

LOT



Tri-Level	72403-1K	2024/07/05
Low	72403-11	2024/07/05
Normal	72403-12	2024/07/05
High	72403-13	2024/07/05

Open vial stability: 14 days ¹

Medonic M51 / Swelab Lumi

Parameter ²	Unit ³	Low		Normal		High	
		Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵
WBC	10 ⁹ /L	3.40	0.50	8.01	1.00	18.20	2.50
Neu%	%	49.3	16.0	54.9	8.0	64.0	8.0
Lym%	%	39.1	21.0	30.0	8.0	20.0	12.0
Mon%	%	6.2	5.0	6.5	5.0	7.0	6.0
Eos%	%	5.4	5.0	8.6	6.0	9.0	8.0
Bas%	%	2.1	10.0	2.0	10.0	1.6	10.0
Neu#	10 ⁹ /L	1.68	0.50	4.40	0.70	11.64	2.00
Lym#	10 ⁹ /L	1.33	0.70	2.40	0.80	3.64	2.00
Mon#	10 ⁹ /L	0.21	0.30	0.52	0.50	1.28	1.10
Eos#	10 ⁹ /L	0.18	0.20	0.69	0.60	1.64	1.40
Bas#	10 ⁹ /L	0.07	0.40	0.16	0.80	0.29	1.80
RBC	10 ¹² /L	2.27	0.24	4.62	0.30	5.41	0.50
HGB	g/L	59	6	135	6	172	8
HCT	%	18.6	2.0	41.6	3.0	52.0	4.0
MCV	fL	81.8	5.0	90.0	5.0	96.2	6.0
MCH	pg	26.0	2.5	29.2	2.5	31.8	2.5
MCHC	g/L	318	30	325	30	331	30
RDW-CV	%	15.1	3.0	15.1	3.0	15.1	3.0
PLT	10 ⁹ /L	51	20	252	40	519	60
MPV	fL	8.8	3.0	8.6	3.0	8.7	3.0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.



1	2	3	4	5
bg	Трайност на отворен флакон: 14 дни	Параметър	Средно	Обхват
cs	Stabilita otevřené ampulky: 14 dní	Parametr	Průměr	Rozsah
da	Åben rør stabilitet: 14 dage	Parameter	Middelværdi	Område
de	Stabilität geöffneter Flaschen: 14 tage	Parameter	Mittelwert	Bereich
el	Διάρθρωση δείγματος μετά την αποσφράγιση: 14 ημέρα	Παράμετρος	Μέση	Αναμενόμενο εύρος α
es	Estabilidad de la cápsula abierta: 14 días	Parámetro	Media	Intervalo
et	Stabiilsus avatud vialis: 14 päevade arv	Parameeter	Keskmine	Vahemik
fr	Stabilité en flacon ouvert: 14 jour	Paramètre	Moyenne	Intervalle
hr	Stabilnost otvorene bočice: 14 dani	Parametar	Sredstvo	Raspon
it	Stabilità della fiala aperta: 14 giorni	Parametro	Media	Intervallo
lt	Stabilumas atidarius buteliuką: 14 dienų	Parametras	Vidurkis	Intervalas
lv	Atvērt flakonu stabilitāti: 14 dienas	Parametrs	Nozīmē	Diapazons
no	Åpen rør/lass stabilitet: 14 dager	Parameter	Middelverdi	Område
pl	Trwałość po otwarciu fiolek: 14 Liczba dni	Kryterium	Srednia	Zakres
pt	Estabilidade após abertura do frasco: 14 dias	Parâmetro	Média	Intervalo
ro	Termenul de valabilitate al fiolei destăcut: 14 zile	Parametru	Medie	Interval
ru	Стабильность вскрытого флакона: 14 дней	Параметр	Среднее	диапазон
sk	Stabilita otvorenej injekčnej liekovky: 14 dní	Parameter	Priemer	Rozsah
sl	Stabilnost odprte stekleničke: 14 dni	Parameter	Povprečno	Območje
sr	Proizvod stabilan nakon otvaranja pakovanja: 14 dana	Parametar	Srednja vrednost	Opseg
sv	Hållbarhet för öppnad flaska: 14 dagar	Enhet	Medelvärde	Intervall
tr	Açık sise stabilitesi gün: 14 gün	Parametre	Ortalama	Aralığı

1504546 31340-12



www.boule.com

Boule
Boule Medical AB
Domnarvsgatan 4
SE-105 53 Spånga, Sweden