

# Certificate of Analysis

**Product** BOULE HYPOCHLORITE 2% CLEANER (500 mL)  
**Product No** 1504113  
**Lot No** 2404-176  
**Manufacturing date** 2024-04-15  
**Expiration Date** 2026-04-14

Test Hypochlorite	Specification
Hypochlorite concentration	2.0 – 2.4 %
pH	11.5 – 12.1
Blank	PLT < 100 x10 <sup>9</sup> /L RBC < 0.05 x 10 <sup>12</sup> /L

Actual values meet specifications.

Date: 2024-04-22



Ulrika Miefalk  
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.  
All information in this document is property of Boule Medical AB.

