

Boule Con-5Diff G2

CONTROL

IVD



| Control | Lot |
|-----------|----------|
| Tri-Level | 62404-1K |
| Low | 62404-11 |
| Normal | 62404-12 |
| High | 62404-13 |

| |
|------------|
| 2024-09-03 |
| 2024-09-03 |
| 2024-09-03 |
| 2024-09-03 |

Open vial stability:

14 days ¹

Quintus

| WBC | 10 ⁹ /L |
|-------|---------------------|
| NEU# | 10 ⁹ /L |
| LYM# | 10 ⁹ /L |
| MON# | 10 ⁹ /L |
| EO# | 10 ⁹ /L |
| BAS# | 10 ⁹ /L |
| NEU | % |
| LYM | % |
| MON | % |
| EO | % |
| BAS | % |
| RBC | 10 ¹² /L |
| Hgb | g/dl |
| Hct | % |
| MCV | fL |
| MCH | pg |
| MCHC | g/dl |
| RDWcv | % |
| Plt | 10 ⁹ /L |
| MPV | fL |

| LOT L | 62404-11 |
|-------------------|----------------------|
| Mean ² | ± Range ³ |
| 4.4 | 0.5 |
| 2.5 | 0.5 |
| 1.5 | 0.4 |
| 0.2 | 0.3 |
| 0.2 | 0.3 |
| 0.0 | 0.2 |
| 56.5 | 9.0 |
| 33.8 | 9.9 |
| 5.5 | 8.7 |
| 4.0 | 4.5 |
| 0.2 | 1.0 |
| 2.28 | 0.15 |
| 5.3 | 0.5 |
| 15.2 | 2.5 |
| 66.6 | 5.5 |
| 23.2 | 2.8 |
| 34.9 | 3.6 |
| 19.1 | 5.0 |
| 83 | 25 |
| 10.6 | 3.0 |

| LOT N | 62404-12 |
|-------------------|----------------------|
| Mean ² | ± Range ³ |
| 8.6 | 0.8 |
| 4.0 | 1.3 |
| 3.7 | 0.8 |
| 0.5 | 0.7 |
| 0.4 | 0.3 |
| 0.0 | 0.2 |
| 46.9 | 8.0 |
| 42.9 | 9.9 |
| 5.7 | 8.5 |
| 4.3 | 3.0 |
| 0.2 | 1.0 |
| 4.21 | 0.18 |
| 11.5 | 0.6 |
| 32.0 | 3.0 |
| 75.9 | 5.5 |
| 27.3 | 2.8 |
| 35.9 | 3.0 |
| 19.4 | 5.0 |
| 246 | 40 |
| 10.5 | 3.0 |

| LOT H | 62404-13 |
|-------------------|----------------------|
| Mean ² | ± Range ³ |
| 21.7 | 2.0 |
| 8.3 | 3.4 |
| 11.3 | 2.7 |
| 0.8 | 1.9 |
| 1.3 | 0.6 |
| 0.0 | 0.2 |
| 38.1 | 8.0 |
| 52.0 | 9.9 |
| 3.8 | 8.3 |
| 5.9 | 3.6 |
| 0.2 | 1.0 |
| 4.73 | 0.30 |
| 17.2 | 0.7 |
| 45.5 | 3.5 |
| 96.1 | 5.5 |
| 36.4 | 2.8 |
| 37.8 | 3.0 |
| 15.7 | 5.0 |
| 525 | 65 |
| 10.8 | 3.0 |

| | CONTROL | N |
|---|----------|----------|
| i | Target X | Target Y |
| | 87 | 112 |

| | |
|----|--|
| bg | Трайност на отворен флакон |
| hr | Stabilnost otvorene bočice |
| cs | Stabilita otevřené ampulky |
| da | Åben rør stabilitet |
| et | Stabiilsus avatud viaalid |
| fr | Stabilité en flacon ouvert |
| de | Stabilität geöffneter Flaschen |
| el | Διατήρηση δείγματος μετά την αποσφράγιση |
| it | Stabilità della fiala aperta |
| lv | Atvērt flakonu stabilitāti |
| lt | Stabilumas atidarius buteliuką |
| no | Åpen rør/glass stabilitet |
| pl | Trwałość po otwarciu fioleki |
| pt | Estabilidade após abertura do frasco |
| ro | Termenul de valabilitate al fiolei desfăcută |
| ru | Стабильность вскрытого флакона |
| sr | Proizvod stabilan nakon otvaranja pakovanja |
| sk | Stabilita otvorenej injekčnej liekovky |
| sl | Stabilnost odprte stekleničke |
| es | Estabilidad de la cápsula abierta |
| sv | Hållbarhet för öppnad flaska |
| tr | Açık şişe stabilitesi |

| | | |
|----------------|------------------|---------------------|
| 14 дни | Средно | Обхват |
| 14 dani | Sredstvo | Raspon |
| 14 dny | Průměr | Rozsah |
| 14 dage | Middelværdi | Område |
| 14 päevade arv | Keskmine | Vahemik |
| 14 jour | Moyenne | Intervalle |
| 14 tage | Mittelwert | Bereich |
| 14 ημέρα | Μέση | Αναμενόμενο εύρος α |
| 14 giorni | Media | Intervallo |
| 14 dienas | Nozīmē | Diapazons |
| 14 d. | Vidurkis | Intervalas |
| 14 dager | Middelverdi | Område |
| 14 Liczba dni | Średnia | Zakres |
| 14 dias | Média | Intervalo |
| 14 zile | Medie | Interval |
| 14 дней | Среднее | диапазон |
| 14 dana | Srednja vrednost | Opseg |
| 14 dní | Priemer | Rozsah |
| 14 dni | Povprečno | Območje |
| 14 dias | Media | Intervalo |
| 14 dagar | Medelvärde | Intervall |
| 14 gün | Ortalama | Aralığı |



www.boule.com

Boule
Boule Medical AB
Domnarvsgatan 4
SE-163 53 Spånga, Sweden

26577-5

1504391

For Quintus with barcode reader

1. In the Main Menu go to "Calibration QC" / "Enter QC reference".
2. Highlight the first Barcode field, scan in each calibrator/control barcode of the column in sequence (1-9) and click "OK".
3. Repeat steps 1 and 2 for each lot to be entered.

Instr.: Quintus
 Lot: 62404-11
 Prod.: Con-diff Low
 ID: 6240411+
 Exp.: 3-Sep-2024

Instr.: Quintus
 Lot: 62404-12
 Prod.: Con-diff Normal
 ID: 6240412+
 Exp.: 3-Sep-2024

Instr.: Quintus
 Lot: 62404-13
 Prod.: Con-diff High
 ID: 6240413+
 Exp.: 3-Sep-2024

