

Certificate of Analysis

Product EXIGO LYSE, RFID (1.9 L)

Product No 1504500 Lot No 2405-190

Manufacturing Date 2024-05-13

Expiration Date 2027-05-12

Test Specification

pН 6.95 - 7.25

195 - 225 Osmolality (mOsm/kg)

Conductivity (mS/cm) 12.10 - 12.80

CBC Performance * Passes test

Particles (PLT/RBC channel) ≤ 5

Actual values meet specifications.

Date:2024-05-17

Pernilla Hammarstedt Rydberg

Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated. All information in this document is property of Boule Medical AB.



^{*} Performance test using calibrated Blood cell counter at Boule Medical reference laboratory. All parameters except mid cells are tested. Accepted reference 95 % to 105 % of expected control blood sample values.