


 Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22404-6K	08/29/24
Low	22404-61	08/29/24
Normal	22404-62	08/29/24
High	22404-63	09/02/24

MEDONIC

Medonic M-series/M32

		CONTROL LOT	L 22404-61
		Mean ²	± Range ³
RBC	10 ¹² /L	2.21	0.12
MCV	fL	67.5	5.0
RDW	%	14.4	4.0
Hct	%	14.9	2.5
Plt	10 ⁹ /L	84	15
MPV	fL	11.6	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	25.3	1.5
MCHC	g/dL	37.6	3.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.6	0.8
Lymph	%	45.0	7.0
Mid	%	5.8	5.0
Gran	%	49.2	7.0

CONTROL LOT	N 22404-62
Mean ²	± Range ³
4.19	0.18
81.3	5.0
14.8	4.0
34.1	3.0
229	30
12.0	1.5
8.2	0.6
11.9	0.4
28.4	1.5
34.9	3.0
3.7	0.6
0.6	0.5
3.9	1.0
44.8	5.0
7.6	5.0
47.6	5.0

CONTROL LOT	H 22404-63
Mean ²	± Range ³
5.11	0.22
92.9	5.0
14.1	4.0
47.5	3.5
517	60
12.0	1.5
20.5	1.8
16.2	0.5
31.7	1.5
34.1	3.0
12.4	1.5
1.2	1.0
6.9	2.5
60.5	5.0
5.7	5.0
33.8	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dny (cs)
 - da Åben rr stabilitet 14 dage (da)
 - et Stabiilsus avatud viaalis 14 pevade (et)
 - fr Stabilit en flacon ouvert 14 jour (fr)
 - de Stabilitt geffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilit della fiala aperta 14 giorni (it)
 - lv Atvrt flakonu stabilitti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuk 14 dienas (lt)
 - no Åpen rr/glass stabilitet 14 dager (no)
 - pl Trwaość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade aps abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cpsula abierta 14 dias (es)
 - sv Hllbarhet fr ppnad flaska 14 dagar (sv)
 - tr Aık ŐiŐe stabilitesi 14 gn (tr)

- 2
- Средно
 - Sredstvo
 - Prmr
 - Middelvrdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μση
 - Media
 - Nozm
 - Vidurkis
 - Middelverdi
 - Srednia
 - Mdia
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvrde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Omrde
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενμενο ερος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Omrde
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralgi



[http://www.medonic.se/
support/downloads](http://www.medonic.se/support/downloads)



1504054

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
 Lot: 22404-61
 Prod.: Con-Diff Low
 ID: 2240461+
 Exp.: 29-Aug-2024

Instr.: M-series
 Lot: 22404-62
 Prod.: Con-Diff Norm
 ID: 2240462+
 Exp.: 29-Aug-2024

Instr.: M-series
 Lot: 22404-63
 Prod.: Con-Diff High
 ID: 2240463+
 Exp.: 2-Sep-2024

<p>1:  Bou l e CBAV2</p>	<p> Bou l e CBAV2</p>	<p> Bou l e CBAV2</p>
<p>2:  2240461+</p>	<p> 2240462+</p>	<p> 2240463+</p>
<p>3:  24082900001153030000</p>	<p> 24082900002153030000</p>	<p> 24090200003153030000</p>
<p>4:  10004030800000000000</p>	<p> 10004030000000000000</p>	<p> 10004030800000000000</p>
<p>5:  22112675500560314925</p>	<p> 41918813501190434130</p>	<p> 51122929501620547535</p>
<p>6:  08415116152531537630</p>	<p> 22930120152841534930</p>	<p> 51760120153171534130</p>
<p>7:  03303015040020301608</p>	<p> 08206037060060503910</p>	<p> 20518124150121006925</p>
<p>8:  14440450700585049270</p>	<p> 14840448500765047650</p>	<p> 14140605500575033850</p>
<p>9:  11119411001905254347</p>	<p> 11110035032536889696</p>	<p> 11115994459265504326</p>