

CERTIFICATE OF ANALYSIS

Product: Boule Cleaning Kit

Kit Lot No: 1986

Product Part #: 501-036

Kit Exp. Date: 2026-03-31

Part Number: 503-085			
Boule Detergent Cleaner Lot # 1983			
Manufacturing Date: 2024-03-14		Quality Testing Date: 2024-03-18	
Parameter	Specification	Result	Pass/Fail
pH (units):	7.60 – 8.40	8.29	Pass

Part Number: 503-086			
Boule Hypochlorite (2%) Cleaner Lot # 1984			
Manufacturing Date: 2024-03-27		Quality Testing Date: 2024-03-28	
Parameter	Specification	Result	Pass/Fail
pH (units):	11.80 – 13.00	12.95	Pass

Part Number: 503-087			
Boule Enzymatic Cleaner Lot # 1985			
Manufacturing Date: 2024-03-25		Quality Testing Date: 2024-04-01	
Parameter	Specification	Result	Pass/Fail
pH (units):	7.60 – 8.40	8.11	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

Allyson
Quality and Regulatory

04/03/24
Date

